

PM4150 Backpack / Carry Bag Instructions

Part #507303

Placing the Portable Oxygen Concentrator (POC) in Carry Bag

- Before placing POC inside carry bag, make sure bag is clean and empty.
- When the POC is properly inserted into carry bag:
 - 1) The POC top will be even with the top of the bag, as shown.
 - 2) The control panel can be seen thru the clear window, as shown.



PRECISION MEDICAL

300 Held Drive
Northampton, PA 18067 USA
Tel: (+001) 610-262-6090
Fax: (+001) 610-262-6080
www.precisionmedical.com
ISO 13485 Certified