

Which POC is right for you?

 PrecisionMedical

Both the Live Active Five® and the EasyPulse POC...

Are powerful, lightweight, and won't get in the way of your life! These portable oxygen concentrators use CMV so patients receive therapeutic oxygen in the first third of inspiration! They both have a 5-year warranty, are FAA approved for use on commercial airlines, and have 5 pulse settings utilizing VPSA technology. Most importantly, both Precision Medical devices consistently produce the highest oxygen purity possible.



live **ACTIVE** Five
Portable Oxygen Concentrator

- More than 6 hours of duration with each battery
- Easy-to-switch external battery
- Bright LCD screen with easy-to-read numbers
- User-replaceable sieve beds
- Custom carry bag with cross-body shoulder strap gives oxygen users a hands-free experience



 **EasyPulse POC**
OXYGEN CONCENTRATOR

- Battery packs with up to 9.4 hours of usage time
- External battery available for extended time
- Fully functional on AC or DC power
- Optional backpack making it comfortable for patients to wear all day long

How do Precision Medical POCs compare to other POCs?

	Oxygen Delivery Settings	Oxygen Concentration*	Weight	Battery Duration	Noise Level at Setting 2	FAA approved	Warranty (on device)
Live Active Five	1-5 Pulse	87%-95.5%	5.0 lbs	Up to 6.5 hours	39.7 dBA	yes	5 years
EasyPulse POC 5	1-5 Pulse	87%-95%	6.6 lbs	Up to 4.7 hours	40 dBA	yes	5 years
Inogen One G5	1-6 Pulse	87%-96%	4.77 lbs with single battery	Up to 4.5 hours with single battery	38 dBA	yes	3 years
SimplyGo Mini	1-5 Pulse	87%-96%	5.0 lbs with standard battery	Up to 4.5 hours with standard battery	43 dBA	yes	3 years
Zen-O Lite	0.5-5 Pulse	87%-96%	5.5 lbs	Up to 4.5 hours	37 dBA	yes	3 years
Caire Freestyle Comfort	1-5 Pulse	87%-95.5%	5.0 lbs with single battery	Up to 8 hours	39.9 dBA	meets FAA guidelines	3 years
AirSep FreeStyle 5	1-5 Pulse	87%-95.5%	6.2 lbs	Up to 5 hours	43 dBA	yes	3 years
iGO2	1-5 Pulse	87%-94%	4.85 lbs	Minimum of 3.5 hours	37.5 dBA	yes	5 years
O2 Concepts Oxlife Liberty	1-9 Pulse 0.5-1.5 Continuous	87%-95%	6.35 lbs	Up to 4.5 hours	40 dBA at 2 l/min	yes	3 years
Belluscura X-PLOR	1-4 Pulse	82%-92%	3.25 lbs with single battery	Up to 4.5 hours with single battery	39 dBA	meets FAA guidelines	3 years

Be cautious of unauthorized online sellers and non FDA approved devices.

*Oxygen purity ranges based on your surroundings. The absolute best that any oxygen unit can achieve is 95.6% in a controlled environment.

