

Removable/External Batteries

The outer packaging requires a UN3480 label.

The label must be printed in color.

Print this label at actual size (4" x 4"). Do not shrink to fit page.

Write a telephone number for additional information on the label.



Non-Removable/Internal Batteries

The outer packaging requires a UN3481 label.

The label must be printed in color.

Print this label at actual size (4" x 4"). Do not shrink to fit page.

Write a telephone number for additional information on the label.



Which Precision Medical product are you returning?



**Live Active Five®
Lithium-ion Battery**
508561



**EasyPulse POC
External Lithium-ion
Battery**
506760



**EasyPulse
Portable Oxygen
Concentrator - 3 Liter**
PM4130P



**EasyPulse
Portable Oxygen
Concentrator - 5 Liter**
PM4150P

Removable/External Batteries

Batteries must be placed in rigid outer packaging and inner packaging designed to prevent damage to the batteries and short circuits (non-conductive materials).

Package must not exceed 66lbs (30kgs).

The outer packaging requires a UN3480 label.

Non-Removable/Internal Batteries

For devices with internal batteries that are not removable, run the battery down to one bar on the battery indicator, then package the device securely.

The outer packaging requires a UN3481 label.

Packaging Guidelines

DO

DO use packaging materials to cushion the battery and prevent it from getting wet.

DO securely pack the battery, so it does not shift inside the shipping box.

DO print the appropriate label in color and adhere it to the shipping box.

DON'T

DON'T ship battery/device without using cushioning, waterproof packaging materials.

DON'T allow the contents to shift inside the shipping box.

DON'T print the appropriate label in black and white.