




The  
**right**  
**connection**  
is **everything.**

**Medical Fittings**

Product Catalog







# Medical Fittings Product Catalog

---

Overview.....	4
Ohmeda Style Quick Connect.....	6
Ohmeda Style Winged Quick Connect.....	7
Compact Ohmeda Style Check Unit (Female Coupler).....	8
Ohmeda Style Check Unit (Female Coupler).....	9
Chemetron Style Quick Connect.....	10
Chemetron Style Check Unit (Female Coupler) .....	11
DISS Female.....	12
DISS Male.....	14
Puritan-Bennett Style Quick Connect .....	16
Puritan-Bennett Style Check Unit (Female Coupler) .....	17
Schrader Style Quick Connect .....	18
Schrader Style Check Unit (Female Coupler).....	19
Medstar Style Quick Connect.....	20
Oxequip Style Quick Connect .....	21
“Y” Block Assemblies and DISS “Y” Assemblies .....	22
“T” Bar Assemblies and DISS “T” Assemblies .....	23
T-Branch, Y-Block, Hex/Quad Block and Rail Blocks.....	26
Extension Manifolds, Pole Clamp and Ferrules.....	27
Medical Grade Hose .....	28
Tubing Nipples .....	29
NPT Pipe Fittings .....	30
Condensation Traps.....	32
Neb Drives and Rescue Flows.....	33
International Fittings.....	34



# The right connection is everything.

## Quality products to make the right connection.


Precision Medical is a global pioneer in respiratory care focused on developing and manufacturing high-quality medical devices. With three main product lines spanning the hospital and homecare markets, our product portfolio covers infants through advanced age and ranges from high-quality oxygen administration to aspiration and suction products.

Precision Medical's durable medical fittings, which are made in the USA, allow full customization and targeting to meet the needs of our customers and users. Our fittings are found in hospitals, surgery centers, and long-term/sub-acute facilities worldwide. We manufacture adapters to fit every need such as both male and female couplers and quick connects for all major styles including Ohmeda, Chemetron, Puritan-Bennett, Schrader, and DISS connections. Customized configurations and hose assemblies are available to meet your unique requirements.

At Precision Medical, we take great pride in providing exceptional customer service. Our knowledgeable representatives are here to help you match the right product to your specific need and have it shipped quickly. We are 100% committed to providing products that are safe, accurate, durable and easy to maintain.

Precision Medical. The Difference is Precision.





# Ohmeda Style Quick Connect



Oxygen



Air



Vacuum



Nitrous Oxide



WAGD

	Oxygen	Air	Vacuum	Nitrous Oxide	WAGD
	2101	2201	2301	2401	2601
1/8" NPT Male					
			2301X		
1/8" NPT Male x 1.5" Long					
	2102	2202	2302	2402	2602
1/8" NPT Female					
	2103	2203	2303	2403	2603
1/4" NPT Male					
	2105				
DISS Male (O2)					
		2205	2305	2405	2605
DISS Male					
	2106	2206	2306	2406	2606
DISS Hex Nut					
	2108	2208	2308	2408	2608
DISS Hand Tight					
	1/4" Hose Barb 2111	1/4" Hose Barb 2211	1/4" Hose Barb 2311 5/16" Hose Barb 2312	1/4" Hose Barb 2411	1/4" Hose Barb 2611 5/16" Hose Barb 2615
Hose Barb					
	1/4" Hose Barb 2111S	1/4" Hose Barb 2211S	1/4" Hose Barb 2311S 5/16" Hose Barb 2312S	1/4" Hose Barb 2411S	1/4" Hose Barb 2611S 5/16" Hose Barb 2615S
Hose Barb Swivel					
	2114				
DISS Male with Check Valve (O2)					
		2214	2314	2414	2614
DISS Male with Check Valve					

# Ohmeda Style Winged Quick Connect

	Oxygen	Air	Vacuum	Nitrous Oxide
	2101W	2201W	2301W	2401W
<b>1/8" NPT Male</b>				
			2301WX	
<b>1/8" NPT Male x 1.5" Long</b>				
	2102W	2202W	2302W	2402W
<b>1/8" NPT Female</b>				
	2103W	2203W	2303W	2403W
<b>1/4" NPT Male</b>				
	2105W			
<b>DISS Male (O2)</b>				
		2205W	2305W	2405W
<b>DISS Male</b>				
	2106W	2206W	2306W	2406W
<b>DISS Hex Nut</b>				
	2108W	2208W	2308W	2408W
<b>DISS Hand Tight</b>				
	<b>1/4" Hose Barb</b> 2111W	<b>1/4" Hose Barb</b> 2211W	<b>1/4" Hose Barb</b> 2311W <b>5/16" Hose Barb</b> 2312W	<b>1/4" Hose Barb</b> 2411W
<b>Hose Barb</b>				
	<b>1/4" Hose Barb</b> 2111WS	<b>1/4" Hose Barb</b> 2211WS	<b>1/4" Hose Barb</b> 2311WS <b>5/16" Hose Barb</b> 2312WS	<b>1/4" Hose Barb</b> 2411WS
<b>Hose Barb Swivel</b>				
	2114W			
<b>DISS Male with Check Valve (O2)</b>				
		2214W	2314W	2414W
<b>DISS Male with Check Valve</b>				



# Compact Ohmeda Style Check Unit

(Female Coupler)



Oxygen



Air



Vacuum



Nitrous Oxide



WAGD



1/8" NPT Female

2160

2260

2360

2460

2660



1/8" NPT Male

2161

2261

2361

2461

2661



1/4" NPT Male

2163

2263

2363

2463

2663



Hose Barb

1/4" Hose Barb  
2164

1/4" Hose Barb  
2264

1/4" Hose Barb  
2364  
5/16" Hose Barb  
2367

1/4" Hose Barb  
2464

1/4" Hose Barb  
2664  
5/16" Hose Barb  
2667



DISS Male (O2)

2165



DISS Male

2265

2365

2465

2665



DISS Hex Nut

2166

2266

2366

2466

2666



DISS Hand Tight

2168

2268

2368

2468

2668



# Ohmeda Style Check Unit

(Female Coupler)

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide	 WAGD
 1/8" NPT Female	2150	2250	2350	2450	2650
 1/8" NPT Male	2151	2251	2351	2451	2651
 1/4" NPT Male	2153	2253	2353	2453	2653
 Hose Barb	1/4" Hose Barb 2154	1/4" Hose Barb 2254	1/4" Hose Barb 2354 5/16" Hose Barb 2357	1/4" Hose Barb 2454	1/4" Hose Barb 2654 5/16" Hose Barb 2657
 DISS Male (O2)	2155				
 DISS Male		2255	2355	2455	2655
 DISS Hex Nut	2156	2256	2356	2456	2656
 DISS Hand Tight	2158	2258	2358	2458	2658

# Chemetron Style Quick Connect



Oxygen



Air



Vacuum



Nitrous Oxide



WAGD



1/8" NPT Male

3101

3201

3301

3401

3601



1/8" NPT Male x 1.5" Long

3301X



1/8" NPT Female

3102

3202

3302

3402

3602



1/4" NPT Male

3103

3203

3303

3403

3603



DISS Male (O2)

3105



DISS Male

3205

3305

3405

3605



DISS Hex Nut

3106

3206

3306

3406

3606



DISS Hand Tight

3108

3208

3308

3408

3608



Hose Barb

1/4" Hose Barb  
3111

1/4" Hose Barb  
3211

1/4" Hose Barb  
3311  
5/16" Hose Barb  
3312

1/4" Hose Barb  
3411

1/4" Hose Barb  
3611  
5/16" Hose Barb  
3615



Hose Barb Swivel

1/4" Hose Barb  
3111S

1/4" Hose Barb  
3211S

1/4" Hose Barb  
3311S  
5/16" Hose Barb  
3312S

1/4" Hose Barb  
3411S

1/4" Hose Barb  
3611S  
5/16" Hose Barb  
3615S



DISS Male with Check Valve (O2)

3114



DISS Male with Check Valve

3214

3314

3414

3614



# Chemetron Style Check Unit

(Female Coupler)



Oxygen



Air



Vacuum



Nitrous Oxide



WAGD



1/8" NPT Female

3150

3250

3350

3450

3650



1/8" NPT Male

3151

3251

3351

3451

3651



1/4" NPT Male


3153

3253

3353

3453

3653



Hose Barb

1/4" Hose Barb  
3154

1/4" Hose Barb  
3254

1/4" Hose Barb  
3354  
5/16" Hose Barb  
3357

1/4" Hose Barb  
3454

1/4" Hose Barb  
3654  
5/16" Hose Barb  
3657



DISS Male (O2)

3155



DISS Male

3255

3355

3455

3655



DISS Hex Nut

3156

3256

3356

3456

3656



DISS Hand Tight

3158

3258

3358

3458







3658

# DISS Female













	Oxygen	Air	Vacuum	Nitrous Oxide	Nitrogen	WAGD	Carbon Dioxide
 Hex Nut- 1/8" NPT with 1" Stem	0106						
 Hex Nut-1/8" NPT with 1" Stem Swivel	0106S						
 Hex Nut- 1/8" NPT with 1 1/2" Stem	0107	0207	0307	0407	0507		504248
 Hex Nut- 1/4" NPT with 1 1/2" Stem					0506	0606	
 Hex Nut- Hose Barb	1/4" Hose Barb 0115	1/4" Hose Barb 0216	1/4" Hose Barb 0316 5/16" Hose Barb 0317	1/4" Hose Barb 0416	1/4" Hose Barb 0516	1/4" Hose Barb 0616 5/16" Hose Barb 0617	1/4" Hose Barb 504154
 Hex Nut- Hose Barb Swivel with 1 1/2" Stem	1/4" Hose Barb 0115S	1/4" Hose Barb 0216S	1/4" Hose Barb 0316S 5/16" Hose Barb 0317S	1/4" Hose Barb 0416S	1/4" Hose Barb 0516S	1/4" Hose Barb 0616S 5/16" Hose Barb 0617S	1/4" Hose Barb 504154S
 Oxygen DISS Female x Oxygen DISS Female	0104						



# DISS Female

	Oxygen	Air	Vacuum	Nitrous Oxide	Nitrogen	WAGD	Carbon Dioxide
 <b>Hand Tight- 1/8" NPT with 1" Stem</b>	0110						
 <b>Hand Tight- Hose Barb Swivel</b>	0110S						
 <b>Hand Tight- 1/8" NPT with 1½" Stem</b>	0111	0209	0309	0409	0509		507883
 <b>Hand Tight- 1/4" NPT with 1½" Stem</b>					0508	0608	
 <b>Hand Tight- Hose Barb</b>	<b>1/4" Hose Barb 0118</b>	<b>1/4" Hose Barb 0218</b>	<b>1/4" Hose Barb 0318 5/16" Hose Barb 0319</b>	<b>1/4" Hose Barb 0418</b>	<b>1/4" Hose Barb 0518</b>	<b>1/4" Hose Barb 0618 5/16" Hose Barb 0619</b>	<b>1/4" Hose Barb 504287</b>
 <b>Hand Tight- Hose Barb Swivel with 1½" Stem</b>	<b>1/4" Hose Barb 0118S</b>	<b>1/4" Hose Barb 0218S</b>	<b>1/4" Hose Barb 0318S 5/16" Hose Barb 0319S</b>	<b>1/4" Hose Barb 0418S</b>	<b>1/4" Hose Barb 0518S</b>	<b>1/4" Hose Barb 0618S 5/16" Hose Barb 0619S</b>	<b>1/4" Hose Barb 504287S</b>

# DISS Male

							
	Oxygen	Air	Vacuum	Nitrous Oxide	Nitrogen	WAGD	Carbon Oxide
<b>DISS Male to 1/8" NPT Male</b>	0701	0731	0751	0771	0791		
<b>DISS Male to 1/4" NPT Male</b>	0703	0733	0753	0773	0792	0794	8946
 <b>1/4" Hose Barb</b>	0709					0601 Demand Check 0602	
 <b>1/4" Hose Barb with Knob</b>	0710	0740	0760	0780	0795		
 <b>1/4" Hose Barb with Knob &amp; Check</b>	0711	0741	0761	0781	0796		
 <b>1/4" Hose Barb with Knob Swivel</b>	0710S	0740S	0760S	0780S	0795S		
 <b>1/4" Hose Barb with Knob &amp; Check Swivel</b>	0711S	0741S	0761S	0781S	0796S		





Oxygen



Air



Vacuum



Nitrous  
Oxide



Nitrogen



WAGD



Carbon Oxide

## Demand Check Valve

DISS Male with  
Demand  
Check Valve  
to 1/8" NPT Male

0712

0742

0782

0797

DISS Male with  
Demand  
Check Valve  
to 1/4" NPT Male

0713

0743

0763

0783

0798

0799

8947

## One Way Check Valve



O2, DISS Male x 1/4" NPT Male

508373



Air, DISS Male x 1/4" NPT Male

508374



1/4" NPT Female x 1/4" NPT Male

508375



Oxygen DISS 90° Elbow  
Male x Hex Nut

1" Stem  
0705  
1½" Stem  
0706



Oxygen Power Take-Off  
(1/8" NPT)

0150



Oxygen DISS Male x  
DISS Male

0704



Oxygen DISS Male x  
1/8" NPT Female

0702



Air Power Take-Off  
(1/8" NPT)

0151

# Puritan-Bennett Style Quick Connect



Oxygen



Air



Vacuum



Nitrous Oxide



WAGD



1/8" NPT Male

4101

4201

4301

4401

4601



1/8" NPT Male x 1.5" Long

4301X



1/8" NPT Female

4102

4202

4302

4402

4602



1/4" NPT Male

4103

4203

4303

4403

4603



DISS Male with Knob (O2)

4105



DISS Male with Knob

4205

4305

4405

4605  
(no knob)



DISS Hex Nut

4106

4206

4306

4406

4606



DISS Hand Tight

4108

4208

4308

4408

4608



Hose Barb with Knob

1/4" Hose Barb  
4111

1/4" Hose Barb  
4211

1/4" Hose Barb  
4311  
5/16" Hose Barb  
4312 (no knob)

1/4" Hose Barb  
4411

1/4" Hose Barb  
4611 (no knob)  
5/16" Hose Barb  
4615 (no knob)



Hose Barb Swivel

1/4" Hose Barb  
4111S

1/4" Hose Barb  
4211S

1/4" Hose Barb  
4311S  
5/16" Hose Barb  
4312S

1/4" Hose Barb  
4411S

1/4" Hose Barb  
4611S  
5/16" Hose Barb  
4615S



DISS Male with Knob and  
Check Valve (O2)

4114



DISS Male with Knob and  
Check Valve

4214

4314

4414

4614  
(no knob)



# Puritan-Bennett Style Check Unit

(Female Coupler)



Oxygen



Air



Vacuum



1/8" NPT Female

4150

4250

4350



1/8" NPT Male

4151

4251

4351



1/4" NPT Male

4153

4253

4353



Hose Barb

1/4" Hose Barb  
4154

1/4" Hose Barb  
4254

1/4" Hose Barb  
4354  
5/16" Hose Barb  
4357



DISS Male (O2)

4155



DISS Male

4255

4355



DISS Hex Nut

4156

4256

4356



DISS Hand Tight

4158

4258

4358

# Schrader Style Quick Connect

	Oxygen	Air	Vacuum	Nitrous Oxide
	6101	6201	6301	6401
<b>1/8" NPT Male</b>				
			6301X	
<b>1/8" NPT Male x 1.5" Long</b>				
	6102	6202	6302	6402
<b>1/8" NPT Female</b>				
	6103	6203	6303	6403
<b>1/4" NPT Male</b>				
	6105			
<b>DISS Male (O2)</b>				
		6205	6305	6405
<b>DISS Male</b>				
	6106	6206	6306	6406
<b>DISS Hex Nut</b>				
	6108	6208	6308	6408
<b>DISS Hand Tight</b>				
	<b>1/4" Hose Barb</b> 6111	<b>1/4" Hose Barb</b> 6211	<b>1/4" Hose Barb</b> 6311	<b>1/4" Hose Barb</b> 6411
<b>Hose Barb</b>				
	6114			
<b>DISS Male and Check Valve (O2)</b>				
		6214	6314	6414
<b>DISS Male and Check Valve</b>				
<b>Schrader Piston Identification</b>				
	<b>Hex Piston (standard)</b>	<b>Swivel Piston (add "S" after part number)</b>		

# Schrader Style Check Unit

(Female Coupler)



Oxygen



Air



Vacuum



Nitrous Oxide



Air/Nitrogen



1/8" NPT Female

6150

6250

6350

6450



1/8" NPT Male

6151

6251

6351

6451



1/4" NPT Male

6153

6253

6353

6453



Hose Barb

1/4" Hose Barb  
6154

1/4" Hose Barb  
6254

1/4" Hose Barb  
6354

1/4" Hose Barb  
6454



DISS Male (O2)

6155



DISS Male

6255

6355

6455



DISS Hex Nut

6156

6256

6356

6456

Schrader Air  
Coupler x  
Nitrogen DISS

6556



DISS Hand Tight

6158

6258

6358

6458

Schrader Air  
Coupler x  
Nitrogen DISS

6558

# Medstar Style Quick Connect



Oxygen



Air



Vacuum



	Oxygen	Air	Vacuum
 <b>1/8" NPT Male</b>	8101	8201	8301
 <b>1/8" NPT Male x 1.5" Long</b>			8301X
 <b>1/8" NPT Female</b>	8102	8202	8302
 <b>1/4" NPT Male</b>	8103	8203	8303
 <b>DISS Male (O2)</b>	8105		
 <b>DISS Male</b>		8205	8305
 <b>DISS Hex Nut</b>	8106	8206	8306
 <b>DISS Hand Tight</b>	8108	8208	8308
 <b>Hose Barb</b>	<b>1/4" Hose Barb</b> 8111	<b>1/4" Hose Barb</b> 8211	<b>1/4" Hose Barb</b> 8311
 <b>Hose Barb Swivel</b>	8111S	8211S	8311S
 <b>DISS Male with Check Valve (O2)</b>	8114MS		
 <b>DISS Male with Check Valve</b>		8214	8314



# Oxequip Style Quick Connect

	Oxygen	Air	Vacuum
 1/8" NPT Male	5101	5201	5301
 1/8" NPT Male x 1.5" Long			5301X
 1/8" NPT Female	5102	5202	5302
 1/4" NPT Male	5103	5203	5303
 DISS Male with Knob	5105		
 DISS Male with Knob		5205	5305
 DISS Hex Nut	5106	5206	5306
 DISS Hand Tight	5108	5208	5308
 Hose Barb	1/4" Hose Barb 5111	1/4" Hose Barb 5211	1/4" Hose Barb 5311
 DISS Male with Knob and Check Valve	5114		
 DISS Male with Knob and Check Valve		5214	5314

# "Y" Block Assemblies and DISS "Y" Assemblies

	Oxygen	Air	Vacuum	Nitrous Oxide	Nitrogen
 <p>Ohmeda Style</p>	2152Y	2252Y	2352Y	2452Y	
with Winged Quick Connect	2152WY	2252WY	2352WY	2352WY	
 <p>Compact Ohmeda Style</p>	2162Y	2262Y	2362Y	2462Y	
with Winged Quick Connect	2162WY	2262WY	2362WY	2462WT	
 <p>Chemetron Style</p>	3152Y	3252Y	3352Y	3452Y	
 <p>Puritan-Bennett Style</p>	4152Y	4252Y	4352Y		
 <p>Schrader Style</p>	6152Y	6252Y	6352Y	6452Y	
 <p>DISS "Y" Assemblies</p>	7211	7212	7213	7214	7240

# “T” Bar Assemblies and DISS “T” Assemblies

	Oxygen	Air	Vacuum	Nitrous Oxide
 <p>Ohmeda Style</p>	2152T	2252T	2352T	2452T
with Winged Quick Connect	2152WT	2252WT	2352WT	2452WT
 <p>Ohmeda Style Offset Right</p>	2152T-R	2252T-R	2352T-R	2452T-R
with Winged Quick Connect	2152WT-R	2252WT-R	2352WT-R	2452WT-R
 <p>Ohmeda Style Offset Left</p>	2152T-L	2252T-L	2352T-L	2452T-L
with Winged Quick Connect	2152WT-L	2252WT-L	2352WT-L	2452WT-L
 <p>Compact Ohmeda Style</p>	2162T	2262T	2362T	2462T
with Winged Quick Connect	2162WT	2262WT	2362WT	2462WT
 <p>Compact Ohmeda Style Offset Right</p>	2162T-R	2262T-R	2362T-R	2462T-R
with Winged Quick Connect	2162WT-R	2262WT-R	2362WT-R	2462WT-R
 <p>Compact Ohmeda Style Offset Left</p>	2162T-L	2262T-L	2362T-L	2462T-L
with Winged Quick Connect	2162WT-L	2262WT-L	2362WT-L	2462WT-L

# “T” Bar Assemblies and DISS “T” Assemblies

	Oxygen	Air	Vacuum	Nitrous Oxide
 <p>Chemetron Style</p>	3152T	3252T	3352T	3452T
 <p>Chemetron Style Offset Right</p>	3152T-R	3252T-R	3352T-R	3452T-R
 <p>Chemetron Style Offset Left</p>	3152T-L	3252T-L	3352T-L	3452T-L
 <p>Puritan-Bennett Style</p>	4152T	4252T	4352T	
 <p>Puritan-Bennett Style Offset Right</p>	4152T-R	4252T-R	4352T-R	
 <p>Puritan-Bennett Style Offset Left</p>	4152T-L	4252T-L	4352T-L	



# “T” Bar Assemblies and DISS “T” Assemblies

	Oxygen	Air	Vacuum	Nitrous Oxide
 <p><b>Schrader Style</b></p>	6152T	6252T	6352T	6452T
 <p><b>Schrader Style Offset Right</b></p>	6152T-R	6252T-R	6352T-R	6452T-R
 <p><b>Schrader Style Offset Left</b></p>	6152T-L	6252T-L	6352T-L	6452T-L
 <p><b>DISS “T” Assemblies</b></p>	7111T	7112T		7114T
 <p><b>DISS “T” Assemblies Offset Right</b></p>	7111T-R	7112T-R		7114T-R
 <p><b>DISS “T” Assemblies Offset Left</b></p>	7111T-L	7112T-L		7114T-L

# T-Branch, Y-Block, Hex/Quad Block and Rail Blocks



**T-Branch**



**Y-Block**



**Hex Block (3 outlet)**



**Quad Block (4 outlet)**

**1/8" NPT Standard**  
7100

**1/8" NPT Offset**  
7100L/R

**1/8" NPT Female**  
7201

**1/4" NPT Female**  
7200

**1/8" NPT Female**  
9054

**1/4" NPT Female**  
9052

**1/4" NPT Female**  
9051

**Oxygen**

**Air**

**Vacuum**



**Hex Block**  
**(1/8" NPT Female Inlet)**

7217

7218

7219



**Quad Block Assembly**  
**(1/8" NPT Female Inlet)**

7220

7221

7222


## **Rail Block Assemblies**

All Rail Block Assemblies include  
Wall Bracket and Plate

7223

7224

7225



**Rail Block Assembly**  
**(DISS Male Inlet)**

## **Replacement Parts**

7226 — Mounting Bracket  
9056 — Mounting Plate  
9054 — Hex Block, 3 outlets  
9058 — Screw (2 required)

*Custom assemblies available upon request.*

## Extension Manifolds, Pole Clamp and Ferrules

## Extension Manifold






## FEATURES:

Economical expansion of single wall outlet.

Permanently mounts to wall.

2 foot conductive hose with 1/8" NPT Female inlet is standard.



Outlets are DISS Male with Check Valves.

	Oxygen	Air	Vacuum
 <p>Extension Manifold-2 Outlet</p>	7230	7233	
 <p>Extension Manifold-2 Outlet with Swivel</p>			7236
 <p>Extension Manifold-3 Outlet</p>	7231	7234	
 <p>Extension Manifold-4 Outlet</p>	7232	7235	
 <p>Oxygen Manifold</p>	505568		

## Pole Clamp

## FEATURES:

1/8" NPT Female holes  
(top and front) - Comes with  
1 plug on bottom

  <b>Pole Clamp</b>	7227A shown	1" Pole 7216
	7227B shown	1-1/2" Pole 7227









# Medical Grade Hose

## FEATURES:









Kink resistant and light weight.  
 FDA medical grade conductive hose.  
 Reinforced with tough nylon braiding.  
 Conveniently packaged in 250 foot rolls or as noted.

Color coded to conform to CGA-C9 standards.  
 Maximum working pressure to 200 psi @ 70°F.  
 Meets NFPA 99 requirements for conductivity (conductive hose only).  
 Keep out of sunlight and in a climate controlled atmosphere

## Non-Conductive Hose

Oxygen	Green		(1/4" ID-250' Roll)	9020-250
Air	Yellow		(1/4" ID-250' Roll)	9021-250
Vacuum	White		(1/4" ID-250' Roll)	9022-250
Nitrous Oxide	Blue		(1/4" ID-250' Roll)	9023-250
Nitrogen	Black		(1/4" ID-250' Roll)	9028-250
Carbon Dioxide	Gray		(1/4" ID-250' Roll)	9033-250
WAGD	Purple		(5/16" ID-250' Roll)	9038-250
	Clear		(1/4" ID-100' Roll)	504202-100

## Conductive Hose

Oxygen	Green		(1/4" ID-250' Roll)	9120-250
Air	Yellow		(1/4" ID-250' Roll)	9121-250
Vacuum	White		(1/4" ID-250' Roll) (5/16" ID-250' Roll)	9122-250 9135-250
Nitrous Oxide	Blue		(1/4" ID-250' Roll)	9123-250
Nitrogen	Black		(1/4" ID-250' Roll)	9128-250
Carbon Dioxide	Gray		(1/4" ID-250' Roll)	9125-250
WAGD	Purple		(1/4" ID-250' Roll)	9124-250
ISO Air	Black & White Striped		(1/4" ID-250' Roll)	9126-250

Custom hose assembly configurator is available at: <http://hoseconfig.precisionmedical.com>



Crimping Tool



Hose Retractor

9130








1/4" Single Clamp	HR-1
1/4" Double Clamp	HR-2
5/16" Single Clamp	HR-3
5/16" Double Clamp	HR-4




# Tubing Nipples


Part Number	Autoclave	Method of Cleaning							
		Ethylene Oxide	Chlorine Dioxide	Hydrogen Peroxide	Gamma Ray	Electron Beam	Vaporized Hydrogen Peroxide (VHP)	Quaternary Cleaners	Sodium Hypochlorite "Bleach Solutions"
0154/0154C		X	X	X	X	X	X	X	X
0158/0158C		X	X	X	X	X	X	X	X
0159/0159C	X	X	X	X	X	X	X	X	X
0159K/0159KC	X	X	X	X	X	X	X	X	X
0115	X	X	X	X	X	X	X	X	X


\* Please Remove O-Ring for Autoclaving \*  
 \* Clean all Tubing Nipples after each patient use \*


 <b>0154 0154C</b> <b>Hand Tight Tubing Nipple</b>	<b>Oxygen</b>	0154
	<b>Oxygen with Chain</b>	0154C
 <b>0158 0158C</b> <b>Wing Nut Tubing Nipple</b>	<b>Oxygen</b>	0158
	<b>Oxygen with Chain</b>	0158C
 <b>0159 0159C</b> <b>Hex Nut Tubing Nipple</b>	<b>Oxygen</b>	0159
	<b>Oxygen with chain</b>	0159C
 <b>0159K 0159KC</b> <b>Knurled Nut Tubing Nipple</b>	<b>Oxygen</b>	0159K
	<b>Oxygen with chain</b>	0159KC
 <b>Disposable Tubing Nipple</b>	<b>Clear</b> <b>Yellow</b> <b>Black</b> <b>White</b> <b>Green</b>	0119 0120 0121 0122 0124
*Note: sold in packs of 25 and 50. Add -25 or -50 at end of part number for those quantities		
 <b>9073 9090</b> <b>Tubing Nipple (for Vacuum Regulator)</b>	<b>Straight</b>	9073
	<b>90°</b>	9090
 <b>Chain Ring Assembly</b>	<b>6" Chain with Closed Ring</b>	8004
	<b>6" Chain with Open Ring</b>	8005


# NPT Pipe Fittings


	1/8" NPT	1/4" NPT
	7030	7031
<b>Pipe Tee</b>		


	7032	7033
<b>Street Tee</b>		

	7034	7035
<b>Female Tee</b>		


	7019	7017
<b>Plug</b>		


	<b>Tee</b> 7081
<b>Drop Ear</b>	<b>Elbow</b> 7082


	7092	
<b>Elbow</b>		


	1/8" NPT	1/4" NPT
	7013	7012
<b>Plug</b>		

	7056	7057
<b>45° Street Elbow</b>		

	7058	7059
<b>90° Elbow</b>		

	7063	7065
<b>90° Street Elbow</b>		

	7064	7045
<b>Master Flow Valve</b>		

	7036	7037
<b>Cross</b>		

# NPT Pipe Fittings

 <p><b>Hex Nipple</b></p>	<p>1/8" NPT Male x 1/8" NPT Male 1/4" NPT Male x 1/4" NPT Male 1/8" NPT Male x 1/4" NPT Male</p>	<p>7020 7021 7022</p>
 <p><b>Connector</b></p>	<p>1/8" NPT Female x 1/8" NPT Male 1/8" NPT Female x 1/4" NPT Male 1/4" NPT Female x 1/8" NPT Male 1/4" NPT Female x 1/4" NPT Male 1/8" NPT Female x 1/8" NPT Male x 2.5" Long</p>	<p>1017 7024 7025 7026 8955</p>
 <p><b>Coupler</b></p>	<p>1/8" NPT Female x 1/8" NPT Female 1/4" NPT Female x 1/4" NPT Female 1/4" NPT Female x 1/8" NPT Female</p>	<p>7027 7028 7029</p>
 <p><b>Hose Barb Connector</b></p>	<p>1/8" NPT Male x 1/8" Diameter 1/8" NPT Male x 3/16" Diameter 1/8" NPT Male x 1/4" Diameter 1/4" NPT Male x 1/4" Diameter 1/8" NPT Male x 5/16" Diameter 1/8" NPT Female x 1/4" Diameter 1/4" NPT Female x 1/4" Diameter 1/4" NPT Female x 5/16" Diameter</p>	<p>7038 7039 7040 7043 8952 7046 7049 8951</p>
 <p><b>Swiveling Hose Barbs</b></p>	<p>1/8" NPT Male x 5/16" Diameter 1/8" NPT Male x 1/4" Diameter</p>	<p>507765 507764</p>
 <p><b>Pipe Nipple</b></p>	<p>1/8" NPT Male x 0.75" Long 1/8" NPT Male x 1.5" Long 1/8" NPT Male x 2" Long 1/8" NPT Male x 2.5" Long 1/4" NPT Male x 1.5" Long 1/4" NPT Male x 2" Long 1/4" NPT Male x 2.5" Long 1/8" NPT Male x 3" Long</p>	<p>7044 7050 7051 7052 7053 7054 7055 7062</p>
 <p><b>Cap</b></p>	<p>1/8" NPT Female</p>	<p>7018</p>
 <p><b>Hose Splicer</b></p>	<p>1/4" Inside Diameter</p>	<p>7041</p>

# Condensation Traps

## Replacement Parts



**Inlet DISS Female  
Oxygen x Outlet DISS  
Male Oxygen**

PM15-4



Jar

507588



**Inlet DISS Female  
Air x Outlet DISS  
Male Air**

PM15-41



Filter

507587



**Inlet DISS Female  
Air x Outlet DISS Male  
Air 90° Down**

PM15-43



Jar and Filter

1271



**Inlet DISS Female  
Oxygen 90° Up x Outlet  
DISS Male Oxygen 90°  
Down**

PM15-44










**Inlet DISS Male  
Air 90° Down x Outlet  
DISS Female Air 90° Up**

PM15-45



# Neb Drives and Rescue Flows

## Neb Drives

		Oxygen	Air
 <p>Ohmeda</p>	6 to 7 l/min	8126	8226
	7 to 8 l/min	8127	8227
	8 to 9 l/min	8128	8228
	10 to 11 l/min	8129	8229
 <p>Chemetron</p>	6 to 7 l/min	8136	8236
	7 to 8 l/min	8137	8237
	8 to 9 l/min	8138	8238
	10 to 11 l/min	8139	8239
 <p>DISS Hand Tight</p>	6 to 7 l/min	8116	8216
	7 to 8 l/min	8117	8217
	8 to 9 l/min	8118	8218
	10 to 11 l/min	8119	8219
 <p>Puritan-Bennett with Knob</p>	6 to 7 l/min	8146	8246
	7 to 8 l/min	8147	8247
	8 to 9 l/min	8148	8248
	10 to 11 l/min	8149	8249
 <p>Schrader</p>	6 to 7 l/min	8166	8266
	7 to 8 l/min	8167	8267
	8 to 9 l/min	8168	8268
	10 to 11 l/min	8169	8269
 <p>Oxequip with Knob</p>	6 to 7 l/min	8156	8256
	7 to 8 l/min	8157	8257
	8 to 9 l/min	8158	8258
	10 to 11 l/min	8159	8259
 <p>Medstar</p>	6 to 7 l/min	8176	8276
	7 to 8 l/min	8177	8277
	8 to 9 l/min	8178	8278
	10 to 11 l/min	8179	8279

## Rescue Flows

	Oxygen 14 to 16 l/min
DISS Hand Tight	8115
Ohmeda	8125
Chemetron	8135
Puritan-Bennett	8145
Oxequip	8155
Schrader	8165
Medstar	8175
	Oxygen 24 to 26 l/min
DISS Hand Tight	8114
Ohmeda	8124
Chemetron	8134
Puritan-Bennett	8144
Oxequip	8154
Schrader	8164
Medstar	8174

# International Fittings

	Oxygen	Air	Vacuum	Nitrous Oxide
FRENCH (AFNOR)				



Probe x 1/8" NPT Male

505979

505980

505981



Probe x 1/4" Hose Barb

505766

505767

505768

## EUROPEAN (NIST)



Nut/Nipple x 1/8" NPT Male

505752

505772

505773



Nut/Nipple x 1/4" Hose Barb

505774

505775

505776



Body Adaptor x 1/8" NPT Male

505777

505778

505779



Body Adaptor x 1/4" Hose Barb

505780

505781

505782

## AUSTRALIAN (SIS)



Nut/Nipple x 1/8" NPT Male

505769

505770

505771

# International Fittings

	Oxygen	Air	Vacuum	Nitrous Oxide
SCANDINAVIAN (AGA)				
 Probe x 1/8" NPT Male	505763	505764	505765	
 Probe x Hose Barb Swivel	1/4" Hose Barb 505763S	1/4" Hose Barb 505764S	1/4" Hose Barb 505765S	
GERMAN (DIN)				
 Probe x 1/8" NPT Male	505753	505761	505762	
 Probe x Hose Barb Swivel	1/4" Hose Barb 505753S	1/4" Hose Barb 505761S	1/4" Hose Barb 505762S	
BRITISH (BOC)				
 Probe x 1/8" NPT Male	6180	6184	6182	6183
 Probe x Hose Barb Swivel	1/4" Hose Barb 6180S	1/4" Hose Barb 6184S	1/4" Hose Barb 6182S	1/4" Hose Barb 6183S
 PK Schrader	6130	6230		



# PrecisionMedical

The Difference is Precision.

300 Held Drive | Northampton, PA 18067

[precisionmedical.com](https://precisionmedical.com)

T: 800-272-7285 | F: 800-352-1240