

*live* **ACTIVE** *Five*®

Portable Oxygen Concentrator



**PrecisionMedical**

| The Difference is **Precision.**

# live **ACTIVE** Five<sup>®</sup>

## Portable Oxygen Concentrator



Lightweight, wearable design weighing only 5 lbs

5 pulse settings utilizing VPSA technology

With CMV, patients receive therapeutic oxygen in the first third of inspiration

Bright LCD screen with large, easy to read numbers

Easily replaceable sieve beds

More than 6 hours of duration with each battery

2 hour battery charge time in the device

Custom carry bag with cross-body shoulder strap gives oxygen users a hands-free experience

Low cost of ownership

## Low cost of ownership

Precision Medical POCs can help your business by eliminating the high recurring costs of deliveries and the burden of managing cylinders. Other manufacturers use basic pressure only systems causing frequent sieve bed replacements. With our dependable POCs you can expect to minimize patient calls and enjoy lower cost of ownership.

## Highest oxygen purity possible

The Live Active Five provides vacuum pressure swing adsorption technology (VPSA) ensuring each sieve bed is rejuvenated every cycle resulting in clean sieve beds and high oxygen purity (95.5%) for a longer period of time.

Controlled Minute Volume (CMV) allows patients to receive therapeutic oxygen in the first third of the inspiratory phase for maximum effectiveness, regardless of breath rate.

## Simple to use with maximum comfort

The device comes with a custom-designed carry bag that allows you to comfortably wear the device on either side of the body or as a cross body. The dual curve design of the device keeps the nasal cannula in place, always keeps the unit's screen and cannula front-facing, and the bag's smart strap prevents it from unhooking from the oxygen outlet.

The plain language on the LCD screen makes the device simple to use. It only takes the push of a button to turn it on and set the desired setting.

## Easy-to-change battery

Each lithium-ion battery packs more than 6 hours of runtime! The top-loading design of the Live Active Five allows you to easily change the battery while the device is still in the carry bag. Additional batteries sold separately adding more security for long trips with the device.



# Ordering Information

---

## Live Active Five Starter Kit

PM4155

includes: Live Active Five POC with battery, AC power adapter with AC power cord, DC car adapter, carry bag, accessory bag and nasal cannula

## Live Active Five Power Bundle

PM4155EB

includes: Live Active Five Starter Kit plus an additional battery

## Live Active Five Battery Charger

508649

includes: battery charger base, AC power adapter with AC power cord and accessory bag

## Live Active Five Battery Charger Bundle

508650

includes: Live Active Five Battery Charger with battery

## Accessories

Live Active Five lithium-ion battery - 508561

Live Active Five AC power adapter with AC power cord - 508660

Live Active Five DC car adapter - 508558

Carry bag with shoulder strap and handle - 508567

Battery charger AC power adapter with AC power cord - 508699

Battery charger DC car adapter - 506750

Sieve bed replacement - 508697



# Product Specifications

---

## Live Active Five

**Oxygen concentration** Up to 95.5%

**Sound level** 39.7 dBA

**Height** 8.4 in (21.4 cm)

**Width** 3.2 in (8.3 cm)

**Depth** 8.5 in (21.6 cm)

**Weight** 5.0 lbs (2.2 kg) with battery

### Alarm types - audible and visual

No breath detected, high breath rate, low battery, blocked cannula, external power failure, low oxygen purity

### User interface

Push buttons with backlit LCD screen

## Battery

**Type** Lithium-ion (94.7 Wh)

**Weight** 1.0 lbs (0.45 kg)

### Battery duration (at any breath rate)

Setting 1 6.5 hrs

Setting 2 4.3 hrs

Setting 3 2.7 hrs

Setting 4 2.0 hrs

Setting 5 1.5 hrs

Specifications are subject to change without prior notice.

