



The Difference is **Precision.**



MEDICAL AIR COMPRESSORS

Medical Air Compressors

Medical Air Compressors are intended to supply a portable compressed air source adequate to drive high volume nebulizers. Precision Medical's reliable Medical Air Compressors deliver hospital-grade performance with safety and maximum durability. The lightweight and compact design ensures easy portability combined with ultra-quiet operation. The long-life induction motor and oil-free pump design eliminates the need for regular maintenance. Used in either a hospital or homecare setting, our durable air compressors deliver the performance and function healthcare providers and caregivers can rely on.

EasyComp Compressor (PM50)

Features

Powerful air flow

Accurate pressure gauge (0-100 psi)

Adjustable regulator knob precisely controls air pressure

Suction cup feet stabilize the unit

Replaceable air filter assembly

Safety built in with high temperature motor shutdown and thermal limit automatic reset

Designed for use in a hospital or homecare setting

Ordering Information

EasyComp Compressor (PM50)

Includes EasyComp Compressor with pressure gauge



EasyComp Compressor is a lightweight, portable medical air compressor that delivers a minimum of 13 l/min of air flow at 50 psi. Weighing only 9.7 lbs, the compact unit has atop facing handle for additional ease. The built-in safety system protects the unit from overheating during operation. Soft, pliable suction cup feet secure the EasyComp to a flat surface and help reduce vibration.

PM50 Accessories:



EasyAir Compressor (PM15)



Features

Continuous operation

Powerful air flow

Pressure control knob for precise control 0-100 psi gauge (PM15-P)

Safety built in with high temperature motor shutdown and thermal limit automatic reset

Designed for use in a hospital or homecare setting

Includes a hospital grade power cord

UL approved

Ordering Information

EasyAir Compressor (PM15-F)

A **flow gauge** sets apart the -F model of the EasyAir Compressor, which includes a J bracket and hospital grade power cord.

EasyAir Compressor (PM15-P)

A **pressure gauge** sets apart the -P model of the EasyAir Compressor, which includes a J bracket and hospital grade power cord.

EasyAir Compressor is a high performance portable medical air compressor designed to reliably supply compressed air 24 hours a day. EasyAir is ideal for continuous use in any clinical or homecare environment. Its sound-dampening case is easy to clean, offers quiet operation and is durable enough for years of heavy use. The oil-free compressor runs cooler and provides clean, breathable air. A precise dial and lock pressure regulator ensures air flow and pressure remain constant. The unit is available with either a pressure gauge or flow gauge.

PM15 Accessories:



Ways to place an order:



Contact your Precision Medical sales representative



Visit precisionmedical.com



Contact Customer Service:
orders@precisionmedical.com
1-800-272-7285

Do you need a special configuration?

Custom configurations, international specifications, special ratings and color options are available.

Contact Customer Service at 1-800-272-7285 for custom requests and more information.